



STEP 2 Push lower edge of ID band down to begin separating imaged band from sheet.

SELF
LAMINATING
WRIST BAND



SELF
LAMINATING
WRIST BAND

STEP 1
Print ID band and
labels on laser printer.

STEP 2
Turn sheet over and lift ID
band from sheet — start from
the right and pull to the left.

STEP 4
Fold film over
imaged area to
create ID band.

THE *St. John*™
COMPANIES

Conf · ID · ent™

Patient Identification System

Reorder: Call (800) 435-4242

Cat. No. WBW23

Licensed under Laser Band U.S. Patent #6,438,881; #6,067,739; #6,000,160; #5,933,993 and others pending and Licensed under U.S. Patent #5,653,472.

