


**STEP 1**  
Print ID band on laser printer

**STEP 2**  
Push down here to begin separating imaged band from sheet.

**STEP 3**  
Turn sheet over and lift ID band from sheet – start from the right and pull to the left.

**STEP 4**  
Fold film over imaged area to create ID band.

SELF LAMINATING WRIST BAND



**WRIST**



THE *A. John* COMPANIES  
**Conf-ID-ent**  
Patient Identification System  
Reorder: Call (800) 435-4242  
Cat. No. WBW24



SELF LAMINATING WRIST BAND

**ANKLE**

Licensed under Laser Band U.S. Patent #6,438,881; #6,067,739; #6,000,160; #5,933,993 and others pending.

Back side is Tamper Evident

