

STEP 2 Push lower edge of ID band down to begin separating imaged band from sheet.

SELF LAMINATING WRIST BAND

SELF LAMINATING WRIST BAND

STEP 1
Print ID band and labels on laser printer.

STEP 2
Turn sheet over and lift ID band from sheet — start from the right and pull to the left.

STEP 4
Fold film over imaged area to create ID band.



Patient Identification System
Reorder: Call (800) 435-4242
Cat. No. WBW26

Licensed under Laser Band® U.S. Patent #4,638,811; #6,067,728; #4,000,100; #5,551,993 and others pending and Licensed under U.S. Patent #5,653,472.

Back side is Tamper Evident

