

## Interested in barcoding your patient ID wristbands?

Bio-Logics patient ID bands can easily be barcoded. Benefit from all of the unique features of the Bio-Logics wristband and easily add a barcode to your patient wristbands at any time!

**There are two different styles of Bio-Logics wristbands that can easily be barcoded:**



1. **Label Based** system which uses a label applied to the back side of the blood bank ID tag. The patient information and barcode(s) are printed directly on a label and the label is applied to the tag (aka card).



2. **Card Based** system uses existing embossers with an imaging module added to print barcodes. We can help you add a barcode to your existing embossed cards.

**Benefits of both Bio-Logics wristband systems above include:**

- A. **24x7 Patient Identification™:** The ability to easily reband your patients immediately at bedside, significantly increasing your compliance percentage while reducing errors associated with unbanded patients.
- B. **Utilize existing printers or embossers:** Label Based wristbands can be laser or thermal printed. Therefore, you may not be required to make any additional investments in new printers and their related ongoing maintenance costs. For Card Based wristbands, most embossers can be easily upgraded with a graphics module to print barcodes if not already equipped.



In-A-Snap™ alert clasps pictured

- C. **Consolidate Multiple Bands:** Because Bio-Logics wristbands integrate the blood band and/or alerts directly into the primary admission band, there is never any confusion or second guessing - all alerts are in one place ... easy to see and replace ... never to be missed. Plus, one band vs. many is MUCH more cost effective! Ask us about our new double-wide, easy to insert Alerts and our brand new In-A-Snap™ alert clasps.



Manual label printer



Electronic printer (part of the ID-Link™ bedside POC system)

- D. **Print Labels On-Demand Directly from the Patient's Wristband:** Prevent ID errors by labeling patient blood samples and/or all patient specimens immediately at bedside with labels printed directly from the patient's wristband. Labels printed on-demand ensures that extra labels aren't floating around and that you don't pay for labels you don't use! Print labels with a manual or an electronic printer at your patient's bedside.
- E. **Bar Code Versatility:** You can print all barcode symbologies, including 2-D, in both vertical and horizontal formatting.
- F. **Improved Barcode Scan Rate:** Larger barcodes printed on the flat ID band surface significantly improve first-scan rate (even on smaller wrists).
- G. **Forward Compatibility:** Rest assured knowing that patient photos, 2-D barcodes and RFID can be easily added at any time.
- H. **Durability of Printed Information:** Designed to withstand the rigors of daily hospital wear with smudge-proof printing and barcoding that is chemical and water resistant even after patient showering.

For more information on barcoding your current patient ID wristband, please contact one of our [patient identification specialists](#).